



## PEDUS-O / PEDUS-U

### Locking Plate System for forefoot osteotomie fixation

#### **PEDUS-O**

**Locking Plate System for fixation of open wedge osteotomies**

With different bar widths, PEDUS-O enables the exact correction of the IM angle with simultaneous stable fixation of the osteosynthesis using locking screws.

#### **PEDUS-U**

**Locking Universal Plate System**

The angular stable universal plate finds its application osteotomy techniques, arthrodesis and fractures in the forefoot and fractures in the forefoot region.

# PEDUS-O

## Open Wedge Locking Plate

# PEDUS-U

## Locking Universal Plate

Fixation of  
Open Wedge Osteotomies



Spacer  
with different levels:  
0, 2, 3, 4, 5, 6 mm


Low Profile  
Protection of the soft tissue environment



Fixation of  
Closing Wedge Osteotomies

Screw holes suitable for:

- Cortical Screws
- Locking Cortical Screws



PEDUS-O / -U Instruments



PEDUS-O WS Open Wedge Plate

12.11121.xxx

Spacer 0, 2, 3, 4, 5, 6 mm

Length 17-23 mm



Locking Cortical Screw Ø 2.3 mm, self-tapping

12.03423.xxx

Length 12-26 mm




PEDUS-U Locking Universal Plate

12.11122.xxx

Hole Distance 12, 14, 16, 18, 20 mm

Length 19-27 mm



Cortical Screw Ø 2.3 mm, self-tapping

12.03424.xxx

Length 12-26 mm

You can also find us on social networks:



Note:  
This description is not sufficient for immediate application of the instruments. Further information is available in the detailed surgical technique.  
Maintenance, care and preparation of Marquardt instruments can be found in the corresponding instructions.