

# TWIN 4.5 & 7.0




## Mid and Rear Foot Screw System

### TWIN

**Two threads, precise placement, stable compression**

Cannulated compression screws in different diameters for a broad range of indications. Compression is achieved by different thread pitches at the screw head and screw tip. The reverse cutting edge facilitates removal.

# TWIN Compression Screws



## TWIN 4.5 / 7.0 Instruments

TWIN 4.5	TWIN 7.0
12.22140.300	12.22140.400



## TWIN Screw Ø 4.5 / 5.5 mm, cannulated

12.03355.xxx(S)
Thread Length 9 - 24 mm
Length 20 - 80 mm



## TWIN Screw Ø 7.0 / 8.5 mm, cannulated

12.03716.xxx(S)	12.03732.xxx(S)
Thread Length 16 mm	Thread Length 32 mm
Length 40 - 120 mm	Length 50 - 120 mm

You can also find us on social networks:



**Note:**

This description is not sufficient for immediate application of the instruments. Further information is available in the detailed surgical technique. Maintenance, care and preparation of Marquardt instruments can be found in the corresponding instructions.